


- Fluorescence/colloidal gold immunochromatographic quantitative analyzer (BXE-QUANT100)



Overall Size: 117 mm × 230 mm × 314 mm

Weight: 1900 g

Major Advantage:

- 1、 High Resolution
- 2、 Good Accuracy
- 3、 Good Repeatability
- 4、 Wide Linear Range
- 5、 Easy To Operate

Normal operating conditions of instrument

- 1、 Ambient temperature range: 10°C-40°C
- 2、 Relative humidity range: $\leq 80\%$
- 3、 Atmospheric pressure range: 86.0kPa ~ 106.0kPa
- 4、 Adapter input 220VAC 50HZ, output 12VDC